

SERIAL NO: _____

St

CHROME MOLY



CE

Customer Contact Details

Customer Name: _____

Advisor: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address:

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Other

VAT Exemption Declaration (UK Only):

I declare that I am chronically sick or disabled by reason of:

and that I am receiving from RGK Wheelchairs Ltd, the wheelchair specified herein and I claim relief from Value Added Tax under group 12 of schedule 8 to Value Added Tax Act 1994.

Signature: _____ Date: _____

RGK Contact Authorisation:

From time to time, RGK would like to send you exclusive information, news and offers on everything RGK, by email or post. We will not pass any of your details to any third party.

I would like to be contacted by RGK:

Signature: _____ Date: _____

SERIAL NO: _____

Guide:

- Base product
- Required option
- Optional Extra
- Critical Information
- Information Needed

SQCXZ010 CLUB SPORT
\$2,089.00

SQCXA148 - BASKETBALL
\$0.00

SQCXA149 - TENNIS *(Only available with XY540 Option 3)*
\$386.00

Max User Weight 125kg (275lbs)

PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM

Standard Specification Includes:

- Removable single sport anti-tip
- Aluminium screw on sideguards
- Adjustable position footrest
- Adjustable centre of gravity (18 degree)
- 72mm castor wheels with extruded forks
- Adjustable height backrest 250-300mm
- Cross laced wire spoke wheels
- Aluminium hard anodised push rim/kenda tyre
- Adjustable seat and backrest upholstery
- 5cm medium density foam cushion (*option 4)
- Calf strap & foot strap

Key Product Measurements

Frame Option (Choose one)

<input type="checkbox"/> XY520 Option 1 <i>(Junior size)</i>	<input type="checkbox"/> XY530 Option 2 <i>(Low Point)</i>	<input type="checkbox"/> XY540 Option 3 <i>(Low/Mid Point)</i>	<input type="checkbox"/> XY550 Option 4 <i>(Max height)</i>
Frame height	Frame height	Frame height	Frame height
<input checked="" type="checkbox"/> Front frame height: 480mm <input checked="" type="checkbox"/> Rear frame height: 420mm	<input checked="" type="checkbox"/> Front frame height: 500mm <input checked="" type="checkbox"/> Rear frame height: 400mm	<input checked="" type="checkbox"/> Front frame height: 500mm <input checked="" type="checkbox"/> Rear frame height: 440mm	<input checked="" type="checkbox"/> Front frame height: 530mm <input checked="" type="checkbox"/> Rear frame height: 530mm
Wheel size	Wheel size	Wheel size	Wheel size
<input type="checkbox"/> 22" XY415 <input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430	<input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430	<input type="checkbox"/> 24" XY420 <input type="checkbox"/> 25" XY430	<input type="checkbox"/> 26" XY440 <input type="checkbox"/> 700c XY450
Seat width (inside of sides)			
<input type="checkbox"/> 340mm <input type="checkbox"/> 370mm	<input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm	<input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm	<input type="checkbox"/> 370mm <input type="checkbox"/> 400mm <input type="checkbox"/> 440mm
Seat length	Seat length	Seat length	Seat length
<input type="checkbox"/> 300mm <input type="checkbox"/> 330mm <input type="checkbox"/> 360mm	<input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm	<input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm	<input type="checkbox"/> 360mm <input type="checkbox"/> 390mm <input type="checkbox"/> 430mm
Frame length	Frame length	Frame length	Frame length
<input checked="" type="checkbox"/> 520mm	<input checked="" type="checkbox"/> 580mm	<input checked="" type="checkbox"/> 580mm	<input checked="" type="checkbox"/> 580mm (SW = 370/400) <input checked="" type="checkbox"/> 630mm (SW = 440)

Options

Finish Options

XN010 Powder coated frame

Specify Code

XN150 Powder coated components (Black)

Accessories

XS020 Sport ratchet strap (click strap) (not fitted) \$81.00

XS010 100mm velcro adjustable strap \$46.00

XK080 Spinerger SLX 24 spoke (Black Spokes) (Wheel upgrade) \$750.00

XP170 Individual backrest embroidery \$136.00
PLEASE COMPLETE PAGE 3

ADDITIONAL NOTES

INDIVIDUAL BACKREST EMBROIDERY

PLEASE PROVIDE JPEG IMAGE AT TIME OF ORDER AND NOTE POSITION BELOW



BRACE BAR

SPEC NOTES: